



# Recombinant Escherichia coli O45:K1 Ribosomal protein S12 methylthiotransferase RimO (rimO)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP486936EOK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>       | B7MGU2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Immunogen Species</b> | Escherichia coli O45:K1 (strain S88 / ExPEC)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Sequence</b>          | MSKVTPQPKI GFVSLGCPKN LVDSEILTE LRTEGYDVVP SYDDADMVIV<br>NTCGFIDSAV QESLEAIGEA LNENGKVIVT GCLGAKEDQI REVHPKVLEI<br>TGPHSYEQVL EHVHHYVPKP KHNPFSLVLP EQGVKLTPRH YAYLKISEGC<br>NHRCTFCIIP SMRGDLVSRP IGEVLSEAKR LVDAGVKEIL VISQDTSAYG<br>VDVKHRTGFH NGEPVKTSMV SLCEQLSKLG IWTRLHYVYP YPHVDDVIPL<br>MAEGKILPYL DIPLQHASPR ILKLMKRP GS VDRQLARIKQ WRKICPELTL<br>RSTFIVGFPG ETEEDFQMLL DFLKEARLDR VGCFKYSPVE GADANALPDQ<br>VPEEVKEERW NRFMQLQQQI SAERLQEKVG REILVIIDEV DEEGAIGRSM<br>ADAPEIDGAV YLNGETNVKP GDILRVKVEH ADEYDLWGSR V |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Target Names</b>      | rimO                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: Ribosomal protein S12 methylthiotransferase RimO Short name= S12 MTTase Short name= S12 methylthiotransferase EC= 2.-.-<br>Alternative name(s): Ribosome maturation factor RimO                                                                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-441                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                            |