



Recombinant Escherichia coli O81 3-methyl-2-oxobutanoate hydroxymethyltransferase (panB)

| | |
|--------------------------|---|
| Product Code | CSB-BP486979EOP |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7MNZ6 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O81 (strain ED1a) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MKPTTIASLQ KCKQDKKRFA TITAYDYSFA KLFADEGINV MLVGDSLGMT VQGHSTLPV TVADIAYHTA AVRRGAPNCL LLADLPFMAY ATPEQAFENA ATVMRAGANM VKIEGGEWLV ETVQMLTERA VPVCGHLGLT PQSVNIFGGY KVQGRGDEAG DRLLSDALAL EAAGAQLLVL ECVPELAKR ITEALAIPI GIGAGNVTDG QILVMHDAFG ITGGHIPKFA KNFLAETGDI RAAVRQYMAE VESGVYPGEE HSFH |
| Source | Baculovirus |
| Target Names | panB |
| Protein Names | Recommended name: 3-methyl-2-oxobutanoate hydroxymethyltransferase EC=2.1.2.11 Alternative name(s): Ketopantoate hydroxymethyltransferase Short name= KPHMT |
| Expression Region | 1-264 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |