



Recombinant Escherichia coli O17:K52:H18 Pyrimidine-specific ribonucleoside hydrolase rihB (rihB)

Product Code	CSB-BP487079EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N5C5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MEKRKIILDC DPGHDDAIAI MMAAKHPVID LLGITIVAGN QTLDKTLING LNVCKLEIN VPVYAGMPQP IMRQQIVADN IHGETGLDGP VFGPLTRQAE STHAVKYIID TLMASDGDIT LVPVGPLSNI AVAMRMQPAI LPKIREIVLM GGAYGTGNFT PSAEFNIFAD PEAARVVFTS GVPLVMMGLD LTNQTVCTPD VIARMERVGG PAGELFSDIM NFTLKTQFEN YGLAGGPVHD ATCIGYLINP DGIKTQEMYV EVDVNSGPCY GRTVCDELGV LGKPANTKVG ITIDTDWFWG LVEECVRGYI KTH
Source	Baculovirus
Target Names	rihB
Protein Names	Recommended name: Pyrimidine-specific ribonucleoside hydrolase rihB EC=3.2.2.8 Alternative name(s): Cytidine/uridine-specific hydrolase
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.