



Recombinant Escherichia coli O17:K52:H18 23S rRNA (uracil (747)-C (5))-methyltransferase RlmC

Product Code	CSB-EP487118EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NAK5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MQCALYDAGR CRSCQWITQP IPEQLSAKTA DLKNLLADFP VEEWCAPVSG PEQGFRNKAK MVSVGSVEKP LLGMLHRDGT PEDLCDCPLY PASFAPVFAA LKPFIARAGL TPYNVARKRG ELKYILLTES QSDGGMMLRF VLRSETKLAQ LRKALPWLQE QLPQLKVITV NIQPVHMAIM EGETEIYLTE QQALADRFND VPLWIRPQSF FQTNPAVASQ LYATARDWVR QLPVNHMWDL FCGVGGFGLH CATPDMQLTG IEIAPEAIAC AKQSAEELGL TRLQFQALDS TQFATAQGEV PELVLVNPFR RGIGKPLCDY LSTMAPRFII YSSCNAQTMA KDIRELPGYR IERVQLFDMF PHTAHYEVLTL LLVKM
Source	E.coli
Target Names	rlmC
Protein Names	Recommended name: 23S rRNA (uracil(747)-C(5))-methyltransferase RlmC EC= 2.1.1.189 Alternative name(s): 23S rRNA(m5U747)-methyltransferase
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.