



Recombinant Escherichia coli O7:K1 4-hydroxy-3-methylbut-2-enyl diphosphate reductase (ispH)

| | |
|--------------------------|---|
| Product Code | CSB-EP487170EON |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7NHD3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O7:K1 (strain IAI39 / ExPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MQILLANPRG FCAGVDRAIS IVENALAIYG APIYVRHEVV HNRVVDLSLR ERGAIFIEQI SEVPDGAILI FSAHGVSQAV RNEAKSRDLT VFDATCPLVT KVHMEVARAS RRGEESILIG HAGHPEVEGT MGQYSNPEGG MYLVESPDDV WKLTVKNEEK LSFMTQTTL S VDDTSDVIDA LRKRFPKIVG PRKDDICYAT TNRQEAVRAL AEQAEVVLVV GSKNSSNSNR LAELAQRMGK RAFLIDDAKD IQEEWVKEVK CVGVTAGASA PDILVQNVVA RLQQLGGGEA IPLEGREENI VFEVPKELRV DIREVD |
| Source | E.coli |
| Target Names | ispH |
| Protein Names | Recommended name: 4-hydroxy-3-methylbut-2-enyl diphosphate reductase EC= 1.17.1.2 |
| Expression Region | 1-316 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |