



Recombinant Escherichia coli O7:K1 Ribosomal RNA small subunit methyltransferase H (rsmH)

Product Code	CSB-BP487173EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NHI8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MMENYKHTTV LLDEAVNGLN IRPDGIYIDG TFGRGGHSRL ILSQFGEEGR LLAIDRDPQA IAVAKTIDDP RFSIIHGPF S ALGDYVAERD LIGKIDGILL DLGVSSPQLD DAERGF S FMR DGPLDMRMDP TRGQSAAEWL QTAAEEADIAW VLKTYGEERF AKRIARAIVE RNREQPMTRT KELAEV VAAA TPVKDKFKHP ATRTFQAVRI WVNSELEEIE QALKSSLNVL APGGRLSIIS FHSLEDRIK RFMRENSRGP QVPAGLPMTE EQLK KLGGRQ LRALGKLMPG EEEVAENPRA RSSVLRIAER TNA
Source	Baculovirus
Target Names	rsmH
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase H EC=2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH
Expression Region	1-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.