



Recombinant Escherichia coli O7:K1 Ribosomal RNA small subunit methyltransferase H (rsmH)

| | |
|--------------------------|---|
| Product Code | CSB-EP487173EON-B |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7NHI8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O7:K1 (strain IAI39 / ExPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MMENYKHTTV LLDEAVNGLN IRPDGIYIDG TFGRGGHSRL ILSQFGEEGR LLAIDRDPQA IAVAKTIDDP RFSIIHGPF S ALGDYVAERD LIGKIDGILL DLGVSSPQLD DAERGF S FMR DGPLDMRMDP TRGQSAAEWL QTAAEEADIAW VLKTYGEERF AKRIARAIVE RNREQPMTRT KELAEV VAAA TPVKDKFKHP ATRTFQAVRI WVNSELEEIE QALKSSLNVL APGGRLSIIS FHSLEDRIK RFMRENSRGP QVPAGLPMTE EQLK KLGGRQ LRALGKLMPG EEEVAENPRA RSSVLRIAER TNA |
| Source | E.coli |
| Target Names | rsmH |
| Protein Names | Recommended name: Ribosomal RNA small subunit methyltransferase H EC=2.1.1.199 Alternative name(s): 16S rRNA m(4)C1402 methyltransferase rRNA (cytosine-N(4)-)-methyltransferase RsmH |
| Expression Region | 1-313 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |