



Recombinant Escherichia coli O7:K1 D-amino acid dehydrogenase small subunit (dadA)

Product Code	CSB-MP487182EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NJF1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRVVILGSGV VGVASAWYLN QAGHEVTVID REPGALET S AANAGQISPG YAAPWAAPGV PLKAIKWMFQ RHAPLAVRLD GTQFQLKMMW QMLRNCDTSH YMENKGRMVR LAEYSRDCLK ALRAETNIQY EGRQGGTLQL FRTEQQYENA TRDIAVLEDA GVPYQLLESS RLAEVEPALA EVAHKLTGGL QLPNDETGDC QLFTQNLARM AEQAGVKFRF NTPVDQLLCD GEQIYGVKCG DEVIKADAYV MAFGSYSTAM LKGIVDIPVY PLKGYSLTIP IAQEDGAPVS TILDETYKIA ITRFDNRIRV GGMAEIVGFN TELLQPRRET LEMVVRDLYP RGGHVEQATF WTGLRPMPD GTPVVGRTF KNLWLNTGHG TLGWTMACGS GQLLSDLLSG RTPAIPYEDL SVARYSRGFT PSRPGHLHGA HS
Source	Mammalian cell
Target Names	dadA
Protein Names	Recommended name: D-amino acid dehydrogenase small subunit EC= 1.4.99.1
Expression Region	1-432
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.