



Recombinant Escherichia coli O7:K1 Probable cytosol aminopeptidase (pepA)

Product Code	CSB-MP487261EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NUH4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MEFSVKSGSP EKQRSACIVV GVFEPRRLSP IAEQLDKISD GYISALLRRG ELEGKPGQTL LLHHVNPVLS ERILLIGCGK ERELDERQYK QVIQKTINTL NDTGSMEAVC FLTELHVKGR NNYWKVRQAV ETAKETLYSF DQLKTNKSEP RRPLRKMVFN VPTRRELTSG ERAIQHGLAI AAGIKAAKDL GNMPPNICNA AYLASQARQL ADSYSKNVIT RVIGEQQMKE LGMHSYLAVG QGSQNESLMS VIEYKGNASE DARPIVLV GK GLTFDSGGIS IKPSEGMDM KYDMCGAAAV YGVMRMVAEL QLPINVIGVL AGCENMPGGR AYRPGDVLTT MSGQTVEVLN TDAEGRLLVC DVLTYVERFE PEAVIDVATL TGACVIALGH HITGLMANHN PLAHELIAAS EQSGDRAWRL PLGDEYQEQL ESNFADMANI GGRPGGAITA GCFLSRFTRK YNWAHLDIAG TAWRSGKAKG ATGRPVALLA QFLLNRAGFN GEE
Source	Mammalian cell
Target Names	pepA
Protein Names	Recommended name: Probable cytosol aminopeptidase EC= 3.4.11.1 Alternative name(s): Leucine aminopeptidase Short name= LAP EC= 3.4.11.10 Leucyl aminopeptidase
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.