



# Recombinant Escherichia coli O127:H6 Probable tRNA threonylcarbamoyladenosine biosynthesis protein Gcp (gcp)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP487299EOB   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7UIX2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O127:H6 (strain E2348/69 / EPEC)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MRVLGIETSC DETGIAIYDD EKGLLANQLY SQVKLHADYG GVVPELASRD<br>HVRKTVPLIQ AALKESGLTA KDIDAVAYTA GPGLVGALLV GATVGRSLAF<br>AWDVPAIPVH HMEGHLLAPM LEDNPPEFPF VALLVSGGHT QLISVTGIGQ<br>YELLGESIDD AAGEAFDKTA KLLGLDYPGG PLLSKMAAQG TAGRFVFRPR<br>MTDRPGLDFS FSGLKTFAAN TIRDNGTDDQ TRADIARAFE DAVVDTLMIK<br>CKRALEQTGF KRLVMAGGVS ANRTLRAKLA EMMKKRRGEV FYARPEFCTD<br>NGAMIAYAGM VRFKAGATAD LGVSVRPRWP LAELPAA        |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | tsaD  |
| <b>Protein Names</b>     | Recommended name: Probable tRNA threonylcarbamoyladenosine biosynthesis protein Gcp<br>Alternative name(s): t(6)A37 threonylcarbamoyladenosine biosynthesis protein   |
| <b>Expression Region</b> | 1-337   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |