



Recombinant Escherichia coli O127:H6 Glucose-1-phosphate adenylyltransferase (glgC)

Product Code	CSB-BP487317EOB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7UKY7
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Purity	>85% (SDS-PAGE)
Sequence	MVSLEKNDHL MLARQLPLKS VALILAGGRG TRLKDLTNKR AKPAVHFGGK FRIIDFALSN CINSGIRRMG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD LLPAQQRMKG ENWYRGTA DA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS HMLIDHVEKG ARCTVACMPV PIEEASAFGV MAVDENDKII EFVEKPANPP SMPNDPSKSL ASMGIVFDA DYLYELLEED DRDENS SHDF GKDLIPKITE AGLAYAHPFP LSCVQSDPDA EPYWRDVGTL EAYWKANLDL ASVPELDMY DRNWPIRTYN ESLPPAKFVQ DRSGSHGRTL NSLVSGGCVI SGSVVVQSVL FSRVRVNSFC NIDSAVLLPE VVVG RSCRLR RCVIDRACVI PEGM VIGENA EEDARRFYRS EEGIVLV TRE MLRKLGHKQE R
Source	Baculovirus
Target Names	glgC
Protein Names	Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27 Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc PPase ADP-glucose synthase
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.