



Recombinant Escherichia coli tRNA modification GTPase MnmE (mnmE)

Product Code	CSB-BP487461ENM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7L851
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain 55989 / EAEC)
Purity	>85% (SDS-PAGE)
Sequence	MSDNDTIVAQ ATPPGRGGVG ILRISGFKAR EVAETVLGKL PKPRYADYLP FKDADGSVLD QGIALWFP GP NSFTGEDVLE LQGHGGPVIL DLLLKRILTI PGLRIARPGE FSERAF LNDK LD LAQAE AIA DLIDASSEQA ARSALNSLQG AFSARVNHLV EALTHLRIYV EAAIDFPDEE IDFLSDGKIE AQLNDVIADL DAVRAEARQG SLLREGMKVV IAGRPNAGKS SLLNALAGRE AAIVTDIAGT TRDVLREHIH IDGMPLHIID TAGLREASDE VERIGIERAW QEIEQADRVL FMVDGTTTDA VDP AEIWPEF IARLPAKLPI TVVRNKADIT GETLGMSEVN GHALIRLSAR TGEGVEVLRN HLKQSMGFD T NMEGGFLARR RHLQALEQAA EHLQQGKAQL LGAWAGELLA EELRLAQQNL SEITGEFTSD DLLGRIFSSF CIGK
Source	Baculovirus
Target Names	mnmE
Protein Names	Recommended name: tRNA modification GTPase MnmE EC= 3.6.-.-
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.