



# Recombinant Escherichia coli Deoxyguanosinetriphosphate triphosphohydrolase (dgt)

<b>Product Code</b>	CSB-EP487518ENM
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7LGM2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain 55989 / EAEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAQIDFRKKI NWHRRYRSPQ GVKTEHEILR IFESDRGRII NSPAIRRLQQ KTQVFPLERN AAVRTRLTHS MEVQQVGRYI AKEILSRLKE LKLEAYGLD ELTGPFESIV EMSCLMHDIG NPPFGHFGEA AINDWFRQRL YPEDAESQPL TDDRCSVAAL RLRDGEPLN ELRRKIRQDL CHFEGNAQGI RLVHTLMRMN LTWAQVGGIL KYTRPAWWRG ETPETHHYLM KKPGYYLSEE AYIARLRKEL NLALYSRFPL TWIMEAADDI SYCVADLEDA VEKRIFTVEQ LYHHLHEAWG QHEKGSLSFSL VVENAWEKSR SNSLSRSTED QFFMYLRVNT LNKLVPYAAQ RFIDNLPAIF AGTFNHALLE DASECSDLLK LYKNVAVKHV FSHPDVEQLE LQGYRVISGL LEIYRPLLSL SLSDFTELVE KERVKRFPPIE SRLFHKLSTR HRLAYVEAVS KLPSDSPEFP LWEYYYRCRL LQDYISGMTD LYAWDEYRRL MAVEQ
<b>Source</b>	E.coli
<b>Target Names</b>	dgt
<b>Protein Names</b>	Recommended name: Deoxyguanosinetriphosphate triphosphohydrolase Short name= dGTP triphosphohydrolase Short name= dGTPase EC= 3.1.5.1
<b>Expression Region</b>	1-505
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.