



Recombinant Escherichia fergusonii Lipoate-protein ligase A (lplA)

Product Code	CSB-MP487562EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LNS6
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MSTLRL LISD SYDPWFNLAV EECIFRQMPA TQRVLFLWRN ADTVVIGRAQ NPWKECNTRR MEEDNVRLAR RSSGGGAVFH DLGNTCFTFM AGKPEYDKTI STSIVLNALN ALGVSAEASG RNDLVVKTAE GDRKVSGSAY RETKDRGFHH GTL LNADLS RLANYLNPKD KKLA AKGITS VRSRVTNLTE LLPGITHEQV CEAITEAFFA HYGERVEAEI ISPDKTPDLP NFAETFARQS SWEWNFGQAP AFSHLLDERF TWGGVELHFD VEKGHITRAQ VFTDSLNPAP LEALAGRLQG CLYRADMLQQ ECEALLVDFP EQEKELRELS AWMAGAVR
Source	Mammalian cell
Target Names	lplA
Protein Names	Recommended name: Lipoate-protein ligase A EC= 2.7.7.63 Alternative name(s): Lipoate--protein ligase
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.