



Recombinant Escherichia fergusonii Sulfate adenylyltransferase subunit 1 (cysN)

Product Code	CSB-MP487618EOR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7LWK4
Product Type	Recombinant Protein
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Purity	>85% (SDS-PAGE)
Sequence	MNTVLAQQIA NEGGVEAWMI AQQHKSLRLF LTCGSVDDGK STLIGRLLHD TRQIYEDQLS SLHNDSKRHG TQGEKLDLAL LVDGLQAERE QGITIDVAYR YFSTEKRKFI IADTPGHEQY TRNMATGAST CDLAILLIDA RKGVLQDQTRR HSFISTLLGI KHLIVAINKM DLVDYSEETF TRIREDYLTAEQLPGNLDI RFVPLSALEG DNVASQSESM LWYSGPTLLE VLETVEIQRV VDAQPMRFPV QYVNRPNLDF RGYAGTLASG RVQVGQRVKV LPSGVESNVA RIVTFDGDRE EAFAGEAITL VLTDEIDISR GDLLAADEV LPAVQSASVA VVWMAEQPLS PGQSYDIKIA GKKTRARVDG IRYQVDINNL TQREVENLSL NGIGLVDLTF DEPLVLDRYQ QNPVTGGLIF IDRLSNVTVG AGMVHEPVSQ ATAAASEFSA FELELNALVR RHFPHWGARD LLGDK
Source	Mammalian cell
Target Names	cysN
Protein Names	Recommended name: Sulfate adenylyltransferase subunit 1 EC= 2.7.7.4 Alternative name(s): ATP-sulfurylase large subunit Sulfate adenylate transferase Short name= SAT
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.