



Recombinant Escherichia coli O8 Pantothenate synthetase (panC)

| | |
|--------------------------|---|
| Product Code | CSB-EP487647E00 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | B7M174 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O8 (strain IA11) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MLIETLPLL RQQIRRLRME GKRVALVPTM GNLHDGHMKL VDEAKARADV VVVSIFVNPM QFDRPEDLAR YPRTLQEDCE KLNKRKVDLV FAPSVKEIYP NGTETHTYVD VPGLSTMLEG ASRPGHFRGV STIVSKLFNL VQPDIACFGE KDFQQLALIR KVMADMGFDI EIVGVPIMRA KDGLALSSRN GYLTAEQRKI APGLYKVLSS IADKLQAGER DLDEIITIAG QELNEKGFRA DDIQIRDADT LLEVSETSKR AVILVAAWLG DARLIDNKMV ELA |
| Source | E.coli |
| Target Names | panC |
| Protein Names | Recommended name: Pantothenate synthetase Short name= PS EC= 6.3.2.1 Alternative name(s): Pantoate--beta-alanine ligase Pantoate-activating enzyme |
| Expression Region | 1-283 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |