



# Recombinant Human N-chimaerin (CHN1)

<b>Product Code</b>	CSB-MP005368HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P15882
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ALTLFDTDE YRPPVWKSYL YQLQQEAPHP RRITCTCEVE NRPKYYGREF HGMISREAAD QLLIVAEGSY LIRESQRQPG TYTLALRFGS QTRNFRLYYD GKHFVGEKRF ESIHDLVTDG LITLYIETKA AEYIAKMTIN PIYEHVGYTT LNREPAYKKH MPVLKETHDE RDSTGQDGVS EKRLTSLVRR ATLKENEQIP KYEKIHNFKV HTFRGPHWCE YCANFMWGLI A QGVKCADCG LNVHKQCSKM VPNDCKPDLK HVKKVYSCDL TTLVKAHTTK RPMVVDMCIR EIESRGLNSE GLYRVSGFSD LIEDVKMAFD RDGEKADISV NMYEDINIIT GALKLYFRDL PIPLITYDAY PKFIESAKIM DPDEQLETLH EALKLLPPAH CETLRYLMAH LKRVTLHEKE NLMNAENLGI VFGPTLMRSP ELDAMAALND IRYQRLVVEL LIKNEIDILF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CHN1
<b>Protein Names</b>	Recommended name: N-chimaerin Alternative name(s): A-chimaerin Alpha-chimerin N-chimerin Short name= NC Rho GTPase-activating protein 2
<b>Expression Region</b>	2-459
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes GTPase-activating protein for p21-rac and a phorbol ester receptor. It plays an important role in ocular motor axon pathfinding. Heterozygous missense mutations in this gene cause Duane s retraction syndrome 2 (DURS2). Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.