



# Recombinant Escherichia coli O8 tRNA sulfurtransferase (thil)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP487669E00-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | B7M3R2  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Escherichia coli O8 (strain IA11)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MKFIILKLFPE ITIKSQSVRL RFIKILTGNI RNVLKHYDET LAVVRHWDNI<br>EVRAKDENQR LAIRDALTRI PGIHHILEVE DVPFTDMHDI FEKALVQYRD<br>QLEGKTFCVR VKRRGKHDFS SIDVERYVGG GLNQHIESAR VKLTNPEVTV<br>HLEVEDDRLL LIKGRYEGIG GFPIGTQEDV LSLISGGFDS GVSSYMLMRR<br>GCRVHYCFFN LGGAAHEIGV RQVAHYLWNR FGSSHRVRFV AINFEPVVGE<br>ILEKIDDGQM GVILKRMMVR AASKVAERYG VQALVTGEAL GQVSSQTLTN<br>LRLIDNVSDT LILRPLISYD KEHIINLARQ IGTEDFARTM PEYCGVISKS<br>PTVKAVKSKI EAEEEEKFDfs ILDKVVEEVN NVDIREIAQQ TEQEVVEVET<br>VNGFGPNdVI LDIRSIDEQE DKPLKVEGID VVSLPFYKLS TKFGDLdQNR<br>TWLLWCERGV MSRLQALYLR EQGFNNVKVY RP |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | thil  |
| <b>Protein Names</b>     | Recommended name: tRNA sulfurtransferase EC= 2.8.1.4 Alternative name(s): Sulfur carrier protein ThiS sulfurtransferase Thiamine biosynthesis protein Thil tRNA 4-thiouridine synthase  |
| <b>Expression Region</b> | 1-482   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.