



Recombinant Escherichia coli O45:K1 Exodeoxyribonuclease 7 large subunit (xseA)

Product Code	CSB-YP487790EOK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MHZ0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O45:K1 (strain S88 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLPSQSPAIF TVSRLNQTVR LLEHEMGQV WISGEISNFT QPASGHWYFT LKDDTAQVRC AMFRNSNRRV TFRPQHGGQV LVRANITLYE PRGDYQIIVE SMQPAGEGLL QLKYEQLKAK LQAEGFLDLQ YKKSLPSPAH CVGVITSKTG AALHDILHVL KRRDPSLPVI IYPTAVQGDD APGQIVRAIE LANQRNECDV LIVGRGGGSL EDLWSFNDER VARAIFASRI PIVSAVGHET DVTIADFVAD LRAPTPSAAA EVVSRNQQEL LRQVQSTHQR LEMAMDYYLA NRTRRFTQIH HRLQQQHPQL RLARQQTMLE RLQKRMSFAL ESQKLRAGQQ QQRLTRQLVQ QNPQSRIHRA QTRIIQLEYR LAETLRAQLS ATRERFGNAV THLEAVSPLS TLARGYSVTS AADGAVLKQV KQVKVGETLT TRLGDGVVIS EVSAVTKTRK SRKKTSNP
Source	Yeast
Target Names	xseA
Protein Names	Recommended name: Exodeoxyribonuclease 7 large subunit EC= 3.1.11.6 Alternative name(s): Exodeoxyribonuclease VII large subunit Short name= Exonuclease VII large subunit
Expression Region	1-458
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.