



# Recombinant Escherichia coli O45:K1 Rhamnulokinase (rhaB)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP487792EOK                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                 |
| <b>Uniprot No.</b>       | B7MI36                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b> | Escherichia coli O45:K1 (strain S88 / ExPEC)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Sequence</b>          | MTFRNCVAVD LGASSGRVML ARYERECRSL TLREIHRFNN GLHSQNGYVT<br>WNVDSLESAL RLGLNKVCEE GIRIDSIGID TWGVDFVLLD QQGQRVGLPV<br>AYRDSRTNGL MAQAQQQLGK RDIYQRSGIQ FLFPNTIYQL RALTEQQPEL<br>IPHIAHALLI PDYFSYRLTG KMNWEYTNAT TTQLVNINSD DWDESLAWS<br>GANKAWFGRP THPGNVIGHW ICPQGNEIPV VAVASHDTAS AVIASPLNGS<br>RAAYLSSGTW SLMGFESQTP FTNDTALAAN ITNEGGAEGR YRVLKNIMGL<br>WLLQRVLQER QINDLPALIA ATQALPACRF IINPNDDRFI NPDEMCSEIQ<br>AACRETAQPI PESDAELARC IFDSLALLYA DVLHELAQLR GEDFSQLHIV<br>GGGCQNTLLN QLCADACGIR VIAGPVEAST LGNIGIQLMT LDELNNVDDF<br>RQVVSTTANL TTFTPNDSE IAHYVAQIHS TRQTKELCA |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | rhaB                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Names</b>     | Recommended name: Rhamnulokinase EC= 2.7.1.5 Alternative name(s):<br>Rhamnulose kinase                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-489                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                        |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                 |