



Recombinant Escherichia coli O81 Formyl-coenzyme A transferase (frc)

Product Code	CSB-BP487877EOP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7MY33
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O81 (strain ED1a)
Purity	>85% (SDS-PAGE)
Sequence	MSTPLQGIKV LDFTGVQSGP SCTQMLAWFG ADVIKIERPG FGDVTRHQLR DIPDIDALYF TMLNSNKRSI ELNTKTAEGK EVMEKLIREA DILVENFHGP AIDHMGFTWE HIQEINPRLI FGSIKGFDEC SPYVNVKAYE NVAQAAGGAA STTGFWDGPP LVSAALGDS NTGMHLLIGL LAALLHREKT GRGQRVTMSM QDAVLNLCRV KLRDQQLDK LGYLEEYPQY PNGTFGDAVP RGGNAGGGGQ PGWILKCKGW ETDPNAYIYF TIQEQNWENT CKAIGKPEWI TDPAYSTAHA RQPHIFDIFA EIEKYTVTID KHEAVAYLTQ FDIPCAPVLS MKEISLDPSL RQSGSVVEVE QPLRGKYLTV GCPMKFSAFT PDIKAAPLLG EHTAAVLQEL GYSDDIEAAM KQNHAH
Source	Baculovirus
Target Names	frc
Protein Names	Recommended name: Formyl-coenzyme A transferase Short name= Formyl-CoA transferase EC= 2.8.3.16
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.