



Recombinant Escherichia coli O17:K52:H18 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase (ispE)

Product Code	CSB-MP487921EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N419
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MRTQWPSPAK LNLFLYITGQ RADGYHTLQT LFQFLDYGDT ISIELRDDGD IRLTPVEGV EHEDNLIVRA ARLLMKTAAAD SGRLPTGSGA DISIDKRLPM GGGLGGGSSN AATVLVALNH LWQCGLSMDE LAEMGLTLGA DVPVFVRGHA AFAEGVGEIL TPVDPPEKWY LVAHPGVSIP TPVIFKDPEL PRNTPKRSIE TLLKCEFSND CEVIARKRFR EVDVAVLSWLL EYAPSRLTGT GACVFAEFD T ESEARQVLEQ APEWLNGFVA KGVNLSPLHR AML
Source	Mammalian cell
Target Names	ispE
Protein Names	Recommended name: 4-diphosphocytidyl-2-C-methyl-D-erythritol kinase Short name= CMK EC= 2.7.1.148 Alternative name(s): 4-(cytidine-5'-diphospho)-2-C-methyl-D-erythritol kinase
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.