



Recombinant Escherichia coli O17:K52:H18 Esterase FrsA (frsA)

Product Code	CSB-BP487953EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N8H0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTQANLSETL FKPRFKHPET STLVRRFNHG AQPPVQSALD GKTIPHWYRM INRLMWIWRG IDPREILDVQ ARIVMSDAER TDDDLYDTVI GYRGGNWIYE WATQAMVWQQ KACAEEDPQL SGRHWLHAAT LYNIAAYPHL KGDDLAEQAQ ALSNRAYEEA AQRLPGTMRQ MEFTVPGGAP ITGFLHMPKG DGPFPVLMC GGLDAMQTDY YSLYERYFAP RGIAMLTIDM PSVGFSSKWK LTQDSSLLHQ HVLKALPNVP WVDHTRVAAF GFRFGANVAV RLAYLES PRL KAVA CLGPVV HTLLSDFKCQ QQVPEMYLDV LASRLGMHDA SDEALRVELN RYS LKVQGLL GRR CPTPMLS GYWKNDPFSP EEDSRLITSS TADGKLEIP FNPVYRNFDK GLQEITDWIE KR LC
Source	Baculovirus
Target Names	frsA
Protein Names	Recommended name: Esterase FrsA EC= 3.1.-.-
Expression Region	1-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.