



Recombinant Escherichia coli O17:K52:H18 Endonuclease 8 (nei)

Product Code	CSB-MP487963EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7N9V3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	PEGPEIRRA ADNLEAAIKG KPLTDVWFAF PQLKSYQSQL IGQHVTHVET RGKALLTHFS NGLTLYSHNQ LYGVWRVVDV GEEPQTTRVL RVKLQTADKT ILLYSASDIE MLTPEQLTTH PFLQRVGPDV LDPNLTPEVV KERLLSPFRF NRQFAGLLLD QAFLAGLGNY LRVEILWQVG LTGNHKAKDL SAAQLDALAH ALLDIPRLSY ATRGQVDENK HHGALFRFKV FHRDGEP CER CGSIIKTTL SSRPFYWCPG CQH
Source	Mammalian cell
Target Names	nei
Protein Names	Recommended name: Endonuclease 8 Alternative name(s): DNA glycosylase/AP lyase Nei EC= 3.2.2.- EC= 4.2.99.18 DNA-(apurinic or apyrimidinic site) lyase Nei Endonuclease VIII
Expression Region	2-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.