



Recombinant Escherichia coli O17:K52:H18 Nicotinate-nucleotide--dimethylbenzimidazole phosphoribosyltransferase (cobT)

Product Code	CSB-BP487978EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NC25
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MQTLADLLNT IP AIDPAAMS RAQRHIDGLL KPVGSLGRLE ALAIQLAGMP GLNGIPHV GK KAVLVMCADH GVWEEGVAIS PKEVTAIQAE NMTRGTTGVC VLAAQAGANV HVIDVGIDTA EPIPGLINMR VARGSGNIAS APAMSRQAE KLLLDVICYT RELAKNGVTL FVGELGMAN TTPAAAIVST ITGRDPEEVV GIGANLPTDK LVNKIDVRR AIMLNQPNPQ DGVDVLAKVG GFDLVGMAGV MLGAASCGLP VLLDGFLSYA AALAACQMSP AIKPYLIPSH LSAEKGARIA LSHLGLEPYL NMEMRLGEGS GAALAMPIIE AACAIYNNMG ELAASNIVLP GNTTSDLNS
Source	Baculovirus
Target Names	cobT
Protein Names	Recommended name: Nicotinate-nucleotide--dimethylbenzimidazole phosphoribosyltransferase Short name= NN:DBI PRT EC= 2.4.2.21 Alternative name(s): N(1)-alpha-phosphoribosyltransferase
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.