



Recombinant Escherichia coli O17:K52:H18 Maltoporin (IamB)

Product Code	CSB-BP488010EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NFY0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	VDFHG YARSGIGWTG SGGEQQCFQT TGAQSKYRLG NECETYAELK LGQEVWKEGD KSFYFDTNVA YSVAQQNDWE ATDPAFREAN VQGKNLIEWL PGSTIWAGKR FYQRHDVHMI DFYYWDISGP GAGLENIDVG FGKLSLAATR SSEAGGSSSF ASNNIYDYTN ETANDVFDVR LAQMEINPGG TLELGVVDYGR ANLRDNYRLV DGASKDGWLF TAEHTQSVLK GFNKFFVQYA TDSMTSQGKG LSQGSGVAFD NEKFAYNINN NGHMLRLDH GAISMGDNWD MMYVGMVQDI NWDNDNGTKW WTVGIRPMYK WTPIMSTVME IGYDNVESQR TGDKNNQYKI TLAQWQAGD SIWSRPAIRV FATYAKWDEK WGYDYNGDSK VNPNYGKAVP ADFNGGSFGR GDSDEWTFGA QMEIWW
Source	Baculovirus
Target Names	IamB
Protein Names	Recommended name: Maltoporin Alternative name(s): Maltose-inducible porin
Expression Region	26-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.