



# Recombinant Escherichia coli O17:K52:H18 Ribosomal RNA small subunit methyltransferase C (rsmC)

<b>Product Code</b>	CSB-MP488016EOG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NH38
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSAFTPASEV LLRHSDDFEQ SRILFAGDLQ DDLPARLDTA ASRAHTQQFH HWQVLSRQMG DNARFSLVAT ADDVADCNTL IYYWPKNKPE AQFQLMNLLS LLPVGTDIFV VGENRSGVRS AEQMLADYAP LNKVDSARRC GLYFGRLEKQ PEFDADKFWG EYSVDGLTVK TLPGVFSRDG LDVGSQLLLL TLTPHTKGKV LDVCGGAGVL SVAFARHSPK IRLTLCDVSA PAVEASRATL AANGVEGEVF ASNVFSEVKG RFDMIISNPP FHDGMQTS LD AAQTLIRGAV RHLNSGGELR IVANAFLPYP DVLDETFGFH EVIAQTGRFK VYRAIMTRQA KKG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	rsmC
<b>Protein Names</b>	Recommended name: Ribosomal RNA small subunit methyltransferase C EC=2.1.1.172 Alternative name(s): 16S rRNA m2G1207 methyltransferase rRNA (guanine-N(2)-)-methyltransferase RsmC
<b>Expression Region</b>	1-343
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.