



Recombinant Escherichia coli O17:K52:H18 Lipoate-protein ligase A (lplA)

Product Code	CSB-EP488017EOG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NH54
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O17:K52:H18 (strain UMN026 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSTLRLISD SYDPWFNLAV EECIFRQMPA TQRVLFLWRN ADTVVIGRAQ NPWKECNTRR MEEDNVRLAR RSSGGGAVFH DLGNTCFTFM AGKPEYDKTI STSIVLNALN ALGVSAEASG RNDLVVKTAE GDRKVSGSAY RETKDRGFHH GTLNADLS RLANYLNPKD KKLAAKGITS VRSRVNTLLE LLPGITHEQV CEAITEAFFA HYGEHVEAEI ISPDKTPDLP NFAETFARQS SWEWNFGQAP AFSHLLDERF TWGGVELHFD VEKGHITRAQ VFTDSLNPAP LEALAGRLQG CLYRADMLQQ ECEALLVDFP EQEKELRELS AWIAGAVR
Source	E.coli
Target Names	lplA
Protein Names	Recommended name: Lipoate-protein ligase A EC= 2.7.7.63 Alternative name(s): Lipoate--protein ligase
Expression Region	1-338
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.