



# Recombinant Escherichia coli O7:K1 3-phenylpropionate-dihydrodiol/cinnamic acid-dihydrodiol dehydrogenase (hcaB)

<b>Product Code</b>	CSB-BP488080EON
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	B7NRJ0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDLHNESIF ITGGGSGGLGL ALVERFIEEG AQVATLELSA AKVASLRQRF GEHILAVEGN VTCYADYQRA LDQILTRSGK LDCFIGNAGI WDHNASLVNT PAETLETGFH ELFNVNVLGY LLSAKACAPA LIASEGSMIF TLSNAAWYYPG GGGPLYTASK HAATGLIRQL AYELAPKVRV NGVGPCGMAS DLRGPQALGQ SETSIMQSLT PEKIAAILPL QFFPQPADFT GPYVMLASRR NNRALSGVMI NADAGLAIRG IRHVAAGLDL
<b>Source</b>	Baculovirus
<b>Target Names</b>	hcaB
<b>Protein Names</b>	Recommended name: 3-phenylpropionate-dihydrodiol/cinnamic acid-dihydrodiol dehydrogenase EC= 1.3.1.87 Alternative name(s): 2,3-dihydroxy-2,3-dihydrophenylpropionate dehydrogenase 3-(cis-5,6-dihydroxycyclohexa-1,3-dien-1-yl)propanoate de
<b>Expression Region</b>	1-270
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.