



Recombinant Escherichia coli O7:K1 L-rhamnose isomerase (rhaA)

Product Code	CSB-MP488103EON
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7NUA3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTTQLEQAW E LAKQRFAAVG IDVEEALRQL DRLPVSMHCW QGDDVSGFEN PEGSLTGGIQ ATGNYPGKAR NASELRADLE QAMRLIPGPK RLNLHAIYLE SDTSVARDQI KPEHFKNWVE WAKANQLGLD FNPSCFSHPL SADGFTLSHA DDSIRQFWID HCKASRRVSA YFGEQLGTPS VMNIWIPDGM KDITVDRLAP RQRLAALDE VISEKLDPAH HIDAVESKLF GIGAESYTVG SNEFYMGYAT SRQTALCLDA GHFHPTEVIS DKISAAMLYV PQLLLHVS RP VRWSDHVVL LDDETQAIAS EIVRHDLFDR VHIGLDFFDA SINRIA AAWVI GTRNMKKALL RALLEPTTEL RKLEAAGDYT ARLALLEEQK SLPWQAVWEM YCQRHDTVPV G SEWLESV RAY EKAILSQRG
Source	Mammalian cell
Target Names	rhaA
Protein Names	Recommended name: L-rhamnose isomerase EC= 5.3.1.14
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.