



Recombinant *Oryza sativa* subsp. *indica* Acyl-[acyl-carrier-protein] desaturase 3, chloroplastic (OsI_07337)

Product Code	CSB-EP488970OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B8AIC3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ASTAGVGG IGNPTPRGKK PFAPWREVPP QVTHTLPPEK KEVFDSLEGW AADTILPYLK PVEESWQPQD HLPDPRSPSF GDEVAALRER AAGLPDDHLV CLVGDMVTEE ALPTYQTMLN TMDGGVRDET GAGGSAWAVW TRAWAAEENR HGDLMNKYLY LTGRVDMRQV EKTIQYLIGS GMDPRTENDP YMGFIYTTFFQ ERATSISHGN TARHAGRHDG AALARVCGTV AADEKRHEAA YAAIVAKLFE VDPDYTVRAF ARMMRRKVAM PARLMYDGAD DRLFARFAAV AQLRGVYTAA DYAGIIEFLV ARWGVPLAA GLSGEGRRAQ DFVCSLGPRF RRMEERAQEA AKRAPPAAAA PFSWIHGRQV QL
Source	<i>E.coli</i>
Target Names	OsI_07337
Protein Names	Recommended name: Acyl-[acyl-carrier-protein] desaturase 3, chloroplastic EC= 1.14.19.2 Alternative name(s): Stearoyl-ACP desaturase 3
Expression Region	33-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.