



Recombinant *Caulobacter vibrioides* Alpha-ketoglutarate-dependent dioxygenase AlkB homolog (alkB)

Product Code	CSB-MP489400DQP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B8GWW6
Product Type	Recombinant Protein
Immunogen Species	<i>Caulobacter crescentus</i> (strain NA1000 / CB15N)
Purity	>85% (SDS-PAGE)
Sequence	MAVVRRAVAA RGLQMIKPL TVVPGFDVWP GLLDISAQRA LVEAVLAGAE QAPFSNYRTA YGKPM SVAMT ALGSLGWTSD ARGYRYVDRH PETGRPW PDM PPALLDLWTV LGDPETPPDS CLVNLYRDGA RMGLHQDRDE ADPRFPVLSI SLGDTAVFRI GGVNRKDPTR SLRLASGDVC RLLGPARLAF HGVDRILPGS SSLVPGGGRI NLTLRRARTA
Source	Mammalian cell
Target Names	alkB
Protein Names	Recommended name: Alpha-ketoglutarate-dependent dioxygenase AlkB homolog EC= 1.14.11.-
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.