



# Recombinant Human Zinc finger protein 7 (ZNF7)

<b>Product Code</b>	CSB-BP026961HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17097
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVVTFGDVA VHFSREEWQC LDPGQRALYR EVMLENHSSV AGLAGFLVFK PELISRLEQG EEPWVLDLQG AEGTEAPRTS KTDSTIRTEN EQACEDMDIL KSESYGTVVR ISPQDFPQNP GFGDVSDSEV WLDSHLGSPG LKVTGFTFQN NCLNEETVVP KTFTKDAQQG CKELGSSGLD CQPLESQGES AEGMSQRCEE CGKGIRATSD IALHWEINTQ KISRCQECQK KLSDCQLQGGH TNNCHGEKPY ECAECGKVFR LCSQLNQHQH IHTGEKPFKC TECGKAFRLS SKLIQHQRH TGEKPYRCEE CGKAFGQSSS LIHHQRIHTG ERPYGCRECG KAFSQQSQLV RHQRTHTGER PYPCKEKGKA FSQSSTLAQH QRMHTGEKAQ ILKASDSPSL VAHQRIHAVE KPFKCDECGK AFRWISRLSQ HQLIHTGEKP YKCNKCTKAF GCSSRLIRHQ RTHHTGEKPFK CDECGKGFVQ GSHLIQHQRH HTGEKPYVCN DCGKAFSQQSS SLIYHQRIHK GEKPYECLQC GKAFSMSTQL TIHQRVHTGE RPYKCNECGK AFSQNSTLFQ HQIIHAGVKP YECSECGKAF SRSSYLIEHQ RIHTRAQWFY EYGNALGEGST FVSRKKVNTI KKLHQCEDCE KIFRWRSHLI IHQRIHTGEK PYKCNDCGKA FNRSSRLTQH QKIHMG
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF7
<b>Protein Names</b>	Recommended name: Zinc finger protein 7 Alternative name(s): Zinc finger protein HF.16 Zinc finger protein KOX4
<b>Expression Region</b>	1-686
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.