



Recombinant Human Zinc finger protein 7 (ZNF7)

Product Code	CSB-EP026961HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17097
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEVVTFGDVA VHFSREEWQC LDPGQRALYR EVMLENHSSV AGLAGFLVFK PELISRLEQG EEPWVLDLQG AEGTEAPRTS KTDSTIRTEN EQACEDMDIL KSESYGTVVR ISPQDFPQNP GFGDVSDSEV WLDSHLGSPG LKVTGFTFQN NCLNEETVVP KTFTKDAQQG CKELGSSGLD CQPLESQGES AEGMSQRCEE CGKGIRATSD IALHWEINTQ KISRQCQCQK KLSDCCLQGKH TNNCHGEKPY ECAECGKVFR LCSQLNQHQH IHTGEKPFKC TECGKAFRLS SKLIQHQRH TGEKPYRCEE CGKAFGQSSS LIHHQRIHTG ERPYGCRECG KAFSQQSQLV RHQRTHTGER PYPCKEKGKA FSQSSTLAQH QRMHTGEKAQ ILKASDSPSL VAHQRIHAVE KPFKCDECGK AFRWISRLSQ HQLIHTGEKP YKCNKCTKAF GCSSRLIRHQ RTHTGEKPFK CDECGKGFVQ GSHLIQHQRH HTGEKPYVCN DCGKAFSQSS SLIYHQRIHK GEKPYECLQC GKAFSMSTQL TIHQRVHTGE RPYKCNECGK AFSQNSTLFQ HQIIHAGVKP YECSECGKAF SRSSYLIEHQ RIHTRAQWFY EYGNALGEGST FVSRKKVNTI KKLHQCEDCE KIFRWRSHLI IHQRIHTGEK PYKCNDCGKA FNRSSRLTQH QKIHMG
Source	E.coli
Target Names	ZNF7
Protein Names	Recommended name: Zinc finger protein 7 Alternative name(s): Zinc finger protein HF.16 Zinc finger protein KOX4
Expression Region	1-686
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.