



# Recombinant Rabbit Tryptophan 5-hydroxylase 1 (TPH1)

<b>Product Code</b>	CSB-BP024100RB
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17290
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Oryctolagus cuniculus (Rabbit)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MIEDNKENKD HSLERGRATL IFSLKNEVGG LIKALKIFQE KHVNLLHIES RKSKRRNSEF EIFVDCDTNR EQLNDIFHLL KSHTNVLSTV PPDNFTMKEE GMESVPWFPK KISDLDHCAN RVLMYGSELD ADHPGFKDNV YRKRRKYFAD LMSYKYGDP IPKVEFTEEE IKTWGTVFRE LNKLYPTHAC REYLKNLPLL SKYCGYREDN IPQLEDISNF LKERTGFSIR PVAGYLSPRD FLSGLAFRVF HCTQYVRHSS DPFYTPEPDT CHELLGHVPL LAEPSFAQFS QEIGLASLGA SEEAQKLAT CYFFTVEFGL CKQDQQLRVF GAGLLSSISE LKHVLSGHAK VKPFDPKITY KQECLITTFQ DVYFVSESEFE DAKEKMREFT KTIKRPFGVK YNPYTRSIQI LKDAKSITNA MNELRHDL DV VSDALGKVSRL QLSV
<b>Source</b>	Baculovirus
<b>Target Names</b>	TPH1
<b>Protein Names</b>	Recommended name: Tryptophan 5-hydroxylase 1 EC= 1.14.16.4 Alternative name(s): Tryptophan 5-monoxygenase 1
<b>Expression Region</b>	1-444
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.