



Recombinant Rabbit Tryptophan 5-hydroxylase 1 (TPH1)

Product Code	CSB-EP024100RB-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17290
Product Type	Recombinant Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥85% (SDS-PAGE)
Sequence	MIEDNKENKD HSLEGRATL IFSLKNEVGG LIKALKIFQE KHVNLLHIES RKSRRNSEF EIFVDCDTNR EQLNDIFHLL KSHTNVLSVT PPDNFTMKEE GMESVPWFPK KISDLHDHCAN RVLMYGSELD ADHPGFKDNV YRKRRKYFAD LMSYKYGDP IPKVEFTEEE IKTWGTVFRE LNKLYPTHAC REYLKNLPLL SKYCGYREDN IPQLEDISNF LKERTGFSIR PVAGYLSPRD FLSGLAFRVF HCTQYVRHSS DPFYTPEPDT CHELLGHVPL LAEPSFAQFS QEIGLASLGA SEEAQKLAT CYFFTVEFGL CKQDGQLRVF GAGLLSSISE LKHVLSGHAK VKPFDPKITY KQECLITTFQ DVYFVSESE DAKEKMREFT KTIKRPFGVK YNPYTRSIQI LKDAKSITNA MNELRHDL DVSDALGKVSRL QLSV
Source	E.coli
Target Names	TPH1
Protein Names	Recommended name: Tryptophan 5-hydroxylase 1 EC= 1.14.16.4 Alternative name(s): Tryptophan 5-monooxygenase 1
Expression Region	1-444
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.