



Recombinant Human C2 calcium-dependent domain-containing protein 4D (C2CD4D)

Product Code	CSB-MP490874HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B7Z1M9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWLLEKAGYK VGAAEPAARW APSGLFSKRR APGPPTSACP NVLTPDRIPQ FFIPRLPDP GGAVPAARRH VAGRGLPATC SLPHLAGREG WAFLPESPHT RRRESLFHGP PPAPAGGLPA AQSRLHVSAP DLRLCRAPDS DTASSPDSSP FGSPRPGLGR RRVSRPHSLS PEKASSADTS PHSPRRAGPP TPPLFHLDLFL CCQLRPTRES VLRLGPRGGQ LRLSTEYQAG PGRLRLRLVS AEGLPRPRSR PGSGGGGCCV VLRLRPRVRP REQQSRVVKC SANPIFNEDF FFDGLGPPDL AARSLRAKVL DRGAGLRRDV LLGECETPLI ALLPPLGGGL GPGSSLAPTH LSL
Source	Mammalian cell
Target Names	C2CD4D
Protein Names	Recommended name: C2 calcium-dependent domain-containing protein 4D
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.