



Recombinant *Oryza sativa* subsp. *indica* Flap endonuclease 1-A (FEN1a)

Product Code	CSB-MP490879OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B8AW67
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MGIKGLTKLL ADNAPKAMKE QKFESYFGRR IAVDASMSIY QFLIVVGRTG METLTNEAGE VTSHLQGMFN RTIRLLEAGI KPVYVFDGKP PDLKKQELAK RYSKREDATK ELTEAVEEGD KDAIEKFSKR TVKVTKQHNE ECKRLLRLMG VPVVEAPCEA EAECAALCIN DMVYAVASED MDSLTFGAPR FLRHLMDPSS KKIPVMEFEV AKVLEEELELT MDQFIDLIL SGCDYCDSEK GIGGQTALKL IRQHGSIESI LENINKDRYQ IPEDWPYQEA RRLFKPNVT LDPELKWNA PDEEGLVQFL VKENGFNQDR VTKAIEKIKF AKNKSSQGRL ESFFKPVVST SVPLKRKDOTS EKPTKAVANK KTKGAGGKKK </p>
Source	Mammalian cell
Target Names	FEN1a
Protein Names	Recommended name: Flap endonuclease 1-A Short name= FEN-1-A EC= 3.1.- .- Alternative name(s): Flap structure-specific endonuclease 1-A
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.