



Recombinant *Oryza sativa* subsp. *indica* Peroxisomal (S)-2-hydroxy-acid oxidase GLO3

Product Code	CSB-EP492165OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B8AU13
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MELITNVSEY EQLAKQKLPK MIYDYYASGA EDQWTLKENR EAFSRILFRP RILIDVSRIN MATNVLGFNI SMPIMIAPSA MQKMAHPEGE LATARAASAA GTIMTLSSWS TSSVEEVNSA APGIRFFQLY VYKDRNIVRQ LVRRRAELAGF KAIALTVDTP RLGRREADIK NRFNLPPHLV LKNFEALDLG KMDKTNDSGL ASYVASQVDR SLSWTDVKWL QTITSLPILV KGVMTAEDTR LAVESGAAGI IVSNHGARQL DYVPATISCL EEVVREAKGR LPVFLDGGVR RGTDFVKALA LGASGVFIGR PVLFSLAVDG EAGVRKVLQM LRDELELTMA LSGCTSLAEI TRNHVITDSD RIRRSRL
Source	E.coli
Target Names	GLO3
Protein Names	Recommended name: Peroxisomal (S)-2-hydroxy-acid oxidase GLO3 EC= 1.1.3.15 Alternative name(s): Glycolate oxidase 3 Short name= GOX 3 Short name= OsGLO3 Short chain alpha-hydroxy acid oxidase GLO3
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.