



Recombinant *Oryza sativa* subsp. *indica* Elongation factor Ts, mitochondrial (EFTS)

Product Code	CSB-MP492814OFF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B8BAI9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	SSEVPASEQM NLIKQLRERT SAPIKDVKAS LVSCNWDIDA AQKDLRKRGV VLAACKSSRT AAEGLLAIQ DEKRAAVVEL NCETDFVARN DVFQYLASSL AKLALSARDP GELVFPFGPD YLEVNLNVNL DHPKLSGETT VQSAVTELA MVGENVKFRF GFIMSTTAHG VVCSYMH TCP QPGLGRLAGL ITLEAEDSNA PLDALQRV GK SIAMHIVATK PLFLSKELVS ASAVENERDI LRTQAESSGK SQMAMEKMVE GRLRKYFEEV VLLEQKYVVN DSTNIKSVLN DLSKEVGSKV TVGNFARMEV GEGVSKLFQQ ILIVTRKRKE WILSK
Source	Mammalian cell
Target Names	EFTS
Protein Names	Recommended name: Elongation factor Ts, mitochondrial Short name= EF-Ts Short name= EF-TsMt
Expression Region	51-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.