



Recombinant Human Uncharacterized protein ENSP00000386791

Product Code	CSB-EP493864HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B9A014
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPRFASPLLR NVIIRSQFDG IKRKQCLQYL KTLRTLQYDG FKTVYFGETN IPESLVTGED ISDGYFIQTP TWCIVHAAGS QGWVPWKYRV FLRDELCIKQ EDSLFSEFCD VVRKAYGKCV IVVKERRQQE EQRPKEDREA EGQFYIPTVI SLASIMCCPE VAKSCGHELL SLPSPCNLYN PLDSAWSSLK WFIINNRRNEF CLQSIDSGYS YQCILFSNLI SKGIERINAS KWRTLTSKVR RWENYYLGKF S
Source	E.coli
Target Names	C21orf140
Protein Names	Recommended name: Uncharacterized protein ENSP00000386791
Expression Region	1-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.