



Alb Antibody, FITC conjugated

Product Code	CSB-PA00010G1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02770
Immunogen	Serum Albumin purified from Rat serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, Alb
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Alb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.