



ALB Antibody, Biotin conjugated

Product Code	CSB-PA00020H1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	I6QYU0
Immunogen	Serum Albumin purified from Pigeon serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Pigeon
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Pigeon
Research Area	Others
Target Names	ALB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.