



Recombinant Arabidopsis thaliana DEAD-box ATP-dependent RNA helicase ISE2, chloroplastic (ISE2), partial

Product Code	CSB-MP495132DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	B9DFG3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Source	Mammalian cell
Target Names	ISE2
Protein Names	Recommended name: DEAD-box ATP-dependent RNA helicase ISE2, chloroplastic EC= 3.6.4.13 Alternative name(s): Protein EMBRYO DEFECTIVE 25 Protein INCREASED SIZE EXCLUSION LIMIT 2 Protein PIGMENT DEFECTIVE 317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.