



# ALB Antibody

<b>Product Code</b>	CSB-PA00030E1Rb
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P19121
<b>Immunogen</b>	Serum Albumin purified from Chicken serum (Native Protein)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Chicken
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.1mol/L NaCl, 15mmol/L NaN <sub>3</sub> , PH7.2
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Serum albumin (Alpha-livetin) (allergen Gal d 5), ALB
<b>Immunogen Species</b>	Chicken
<b>Research Area</b>	Others
<b>Target Names</b>	ALB
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.