



# ALB Antibody

Product Code	CSB-PA00030E1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19121
Immunogen	Serum Albumin purified from Chicken serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.1mol/L NaCl, 15mmol/L NaN <sub>3</sub> , PH7.2
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (Alpha-livetin) (allergen Gal d 5), ALB
Immunogen Species	Chicken
Research Area	Others
Target Names	ALB