



ALB Antibody, Biotin conjugated

Product Code	CSB-PA00030H1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19121
Immunogen	Serum Albumin purified from Chicken serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Chicken
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin (Alpha-livetin) (allergen Gal d 5), ALB
Immunogen Species	Chicken
Research Area	Others
Target Names	ALB