



ALB Antibody, HRP conjugated

Product Code	CSB-PA00060F1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02768
Immunogen	Serum Albumin purified from Human serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ALB