



ALB Antibody

Product Code	CSB-PA00060E1Rb
Abbreviation	Serum Albumin purified from serum
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02768
Immunogen	Serum Albumin purified from Human serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Serum albumin, the main protein of plasma, has a good binding capacity for water, Ca(2+), Na(+), K(+), fatty acids, hormones, bilirubin and drugs. Its main function is the regulation of the colloidal osmotic pressure of blood. Major zinc transporter in plasma, typically binds about 80% of all plasma zinc.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, ALB
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ALB