



Alb Antibody, Biotin conjugated

Product Code	CSB-PA00080H1Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07724
Immunogen	Serum Albumin purified from Mouse serum (Native Protein)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serum albumin, Alb, Alb-1 Alb1
Immunogen Species	Mus musculus (Mouse)
Research Area	Cardiovascular
Target Names	Alb
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.